

Ordering Information



- Guidewire compatibility: **0.035"**
- Nominal Pressure: **8atm** (All Models)
- Rated Burst Pressure: **14atm** (All Models)
- Catheter length: **130cm** (All Models)

Balloon Diameter	Balloon Length			No. of Scoring Elements	Compatible Introducer Sheath
	40mm	60mm	80mm		
8.0mm	DLP-8040	DLP-8060	DLP-8080	5	9Fr
9.0mm	DLP-9040	DLP-9060	DLP-9080	5	9Fr
10.0mm	DLP-10040	DLP-10060	DLP-10080	5	9Fr
12.0mm	DLP-12040	DLP-12060	DLP-12080	6	10Fr
14.0mm	DLP-14040	DLP-14060	DLP-14080	6	10Fr

D•Kutting™ XL



PTA Scoring Balloon Dilatation Catheter XL

CE MDR: 816124



DK Medical Technology Co.,Ltd.
☎ +86-512-81880801
✉ cs@dkmedtech.com
📍 Unit #5 Xinhong Park, No. 9 West Suhong Road, Suzhou



Precise Dilatation, Minimized Complication

D•Kutting™ XL

PTA scoring balloon dilatation catheter XL

Go HIGH

Non-compliant balloon with RBP up to 14atm, delivering controlled force where it counts

Go with CONTROL

Triangular scoring elements, fixed and 0.36mm high, for targeted, predictable plaque modification

Go BIG

Largest scoring balloon available, with diameters up to 14mm

