



NEUROSAFE

Glutton™
Aspiration Catheter




Large Lumen
Ultra Flexibility
Great Efficiency


NeuroSafe Medical Co., Ltd.

Glutton™ Aspiration Catheter

The Glutton™ Aspiration Catheter is intended for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (with the internal carotid, middle cerebral – M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset.



Hybrid Structure of Coil and Braiding,
BALANCED KINK RESISTANCE AND TRACKABILITY.



Larger inner diameter **0.073"**

Larger lumen, Stronger TRF, Higher FPE

FEATURES & BENEFITS

- - Hybrid structure of coil and braiding, combined flexibility and maneuverability;
- - Larger inner diameter lumen, higher aspirating efficiency;
- - Distal multi-segment design of gradual hardness;
- - Proximal composite polymer enhanced pushability;
- - Compatible with 6F long sheath;
- - Hydrophilic coating on the distal segment easily reach distal vessels / occlusions.

Order Information:

Catalog No.	Fr. Size	Inner Diameter (inch)	Outer Diameter (inch)	Effective Length (cm)
5F-115	5F	0.058	0.070	115
5F-125	5F	0.058	0.070	125
5F-132	5F	0.058	0.070	132
6F-105	6F	0.073	0.085	105
6F-115	6F	0.073	0.085	115
6F-125	6F	0.073	0.085	125
6F-115BGC	6F	0.070	0.082	115
6F-125BGC	6F	0.070	0.082	125
6F-132BGC	6F	0.070	0.082	132
6F-115S	6F	0.073	0.085	115
6F-125S	6F	0.073	0.085	125
6F-132S	6F	0.073	0.085	132

NeuroSafe Medical Co., Ltd.

Address: Building B, No. 10, Keji 1st Rd, Hi-Tech Industrial Development Zone, Zhuhai, Guangdong, China

Tel: +86 (756) 3635 345 | Email: sales@neurosafe.cn | Postcode: 519 000

Website: www.neurosafe.cn