




Rotational Atherectomy & Thrombectomy Device





Unique catheter tip configuration
for removal of calcified material

No need to pass the occlusion
with a guidewire

Fast and easy set up

No capital equipment required



WATCH VIDEO



ROTATIONAL ATHERECTOMY & THROMBECTOMY DEVICE

BYCROSS is a 6F bi-directional rotational atherectomy & thrombectomy device that is inserted into the peripheral arterial vascular system via a guidewire and an introducer sheath to restore chronically occluded vessels.

- ● ● The BYCROSS catheter consists of a coaxial, flexible, rotating shaft with an expandable tip. The BYCROSS device is battery powered. The calcified material is aspirated via an integrated pump and discharged into a collection bag.



BYCROSS is a trademark of TARYAG Medical Ltd.

Germany

plus medica GmbH & Co. KG
Willstätterstraße 13
40549 Düsseldorf
Germany
Tel: +49 211 54 59 08 80
kundenservice@plusmedica.de
info@plusmedica.de
www.plusmedica.de

Austria

plus medica GmbH
Klingerstraße 9
2353 Guntramsdorf
Austria
Tel: +43 1 50 32 585
kundenservice@plusmedica.at
info@plusmedica.at
www.plusmedica.at

Switzerland

plus medica AG
Oberwilerstraße 3
8442 Hettlingen
Switzerland
Tel: +41 52 316 44 36
kundenservice@plusmedica.ch
info@plusmedica.ch
www.plusmedica.ch

Sweden

plus medica Nordic AB
Lövsjögsvägen 12
70376 Örebro
Sweden
Tel: +46 70 588 70 91
gzivkovic@plusmedica.se
www.plusmedica.se



TARYAG Medical Ltd.

14 Ha'Ilan St. Or Akiva,
3065101 Israel
Tel: +972.4.6776929
Fax: +972.77.3180728
info@taryag-medical.com
www.taryag-medical.com

FACTS & FIGURES

Art. No.	Description	VE
BC213-50	BYCROSS 6F 50CM	1
BC213-70	BYCROSS 6F 70CM	1
BC213-95	BYCROSS 6F 95CM	1

In combination with other devices / not included:

Device	Description
Guidewire	Standard 0,035" (0,89 mm) Guidewire in appropriate length
Syringe	Standard-Luer-Syringe
Guiding Sheath	Stainless Steel reinforced 6F / 7F Guiding Sheath in appropriate length

